OPERATING INSTRUCTIONS

Humidity Cabinet

39-1300, 39-1601

<p>ELE International Chartmoor Road, Chartwell Business Park Leighton Buzzard, Bedfordshire, LU7 4WG England phone: +44 (0) 1525 249200 fax: +44 (0) 1525 249249 email: ele@eleint.co.uk http://www.ele.com ELE International, a division of Hach Lange Ltd.</p>	<p>Distributor:</p>	<p>ELE International Soiltest Product Division PO Box 389, Loveland, CO 80539 USA phone: +1 (800) 323 1242 fax: +1 (970) 663 9781 email: soiltest@eleusa.com http://www.eleusa.com</p>
<p><i>In the interests of improving and updating its equipment, ELE reserves the right to alter specifications to equipment at any time.</i> ELE International 2018 ©</p>		

Contents

	Section	Page
1	General Description	3
2	Installation	4
3	Temperature Setting Procedure	4
4	Humidity Control	7
5	Maintenance	7
	APPENDIX A : Digital Data Recorder (Temperature & Humidity)	9
	APPENDIX B : Programming an Offset Parameter into a JUMO Dtron 316 temperature controller	18
	APPENDIX C : TLK38 OFFSET Instructions	19
	WEEE Directive	

1 General Description

The cooled laboratory humidity cabinet provides the user with extremely accurate and safe temperature control over a wide temperature range.

1.1 Heating

Low wattage heaters are attached to the outer surface of the inner chamber in such a manner as to ensure even heating throughout the chamber.

1.2 Cooling

Direct expansion coils are coupled to a hermetic condensing unit via a solenoid valve. When the refrigeration switch is switched on, the compressor will run continuously. The coolant flow is controlled by a solenoid valve, which is electronically opened and closed by the temperature controller.

1.3 Temperature Controller

The temperature is accurately controlled by a Jumo dTRON 316 microprocessor temperature controller. The controller uses the signal from a PT100 sensor located in the chamber to control the current supplied to the heaters. The temperature can be easily adjusted. See Section 3 for details.

1.4 Over / Under Temperature Safety Cut Out

A TLK38 digital cut out protects the contents of your incubator by allowing you to program over & under temperature safety cut out set points. When the chamber temperature goes out of limits, an audible alarm and visible red neon will bring to your attention any problems. Reset buttons are provided below the cut out for you to manually reset the out-of-limit system.

1.5 Circulating Fan

Air is drawn up through the chamber by a circulating fan located in the top and returned to the base via a full width duct at the rear of the chamber. A door micro-switch stops the fan when the door is opened.

1.6 Humidity

A wide range of fixed relative humidities can be maintained by adding a suitable saturated salt solution to the base of the chamber. See Section 4 for details.

2 Installation

2.1 Connect the 3-core cable to a 240V AC, 50Hz mains supply as follows:

BROWN	-	LIVE
BLUE	-	NEUTRAL
GREEN/YELLOW	-	EARTH

2.2 The incubator is internally protected by a fuse, which isolates the cabinet from the electrical supply if excessive current is drawn. The fuse holder is located on the rear of the cabinet. The instruments on the control panel will illuminate when the supply is connected correctly and POWER switch is on.

2.3 General Precautions

A qualified electrician or other competent person must carry out any electrical work required to install the incubator.

Do not place the incubator in direct sunlight or near a heat source.

Make sure the incubator is not standing on its electrical supply cable.

Before the incubator is switched on for the first time, it should stand still for 2 hours after final positioning to allow the refrigerant to settle.

Keep all ventilation grilles clear. Failure to do so may result in component failure.

The refrigeration system gives off heat when operating. Always ensure there is adequate ventilation around the incubator.
NEVER PLACE ANYTHING ON TOP OF THE INCUBATOR.

Before any cleaning or maintenance work is carried out, the mains supply must be switched off and the plug removed from the electrical socket.

3 Temperature Setting Procedure

Switch both the POWER and COOLING switches ON. DO NOT switch the HIGH HUMIDITY switch on yet! *Please see Section 4.1 for more details.*

3.1 The JUMO dTRON 316 microprocessor temperature controller on your incubator has two LED displays. The upper display, which is red, represents the actual chamber temperature. The lower display, which is green, represents the target temperature. The controller has 4 buttons marked as follows;

- PGM** - Program button (Enters the programming mode).
- ▲** - Increases the temperature set point
- ▼** - Decreases the temperature set point.
- EXIT** - Exits the programming mode.

The **yellow 3** on the temperature controller shows when heat is being supplied. The **yellow 5** shows when cooling is being supplied.

3.2 The incubator operating temperature or set point (SP) can be altered by simply pressing either the **▲** button to increase the temperature or the **▼** button to decrease the temperature. When the required temperature has been entered, release all buttons and the green display will blink once to confirm that a new target temperature has been stored.

3.3 Over / Under Temperature Safety Cut Out

A TLK38 digital safety cut out unit provides essential over & under temperature protection for the products inside the chamber. The cut out has buttons labelled as follows:

- P** - Enters the programming mode.
- U** - *No function.*
- ▲** - **INCREASES** the displayed value.
- ▼** - **DECREASES** the displayed value.

3.4 Setting An Over Temperature Cut-Out Value

If, for example, you would like to operate your incubator at +20.0°C, an over temperature cut-out value needs to be programmed into the cut-out. We recommend that the cut-out value is set to **2 degrees higher** than your chamber temperature. In the example above, the over temperature cut-out set point needs to be set +22°C. To do this follow steps 1 - 3 below:

- a) Press the **[P]** button **once only**. Flashing **AL1** will be displayed plus a value (the **over temperature** cut out value).
- b) Press the **[▲]** button to increase the **over temperature** cut-out value or the **[▼]** button to decrease the **over temperature** cut-out value.
- c) Once the value of your choice is displayed, press the **[P]** button again and the new value will be stored and retained in permanent memory.

3.5 Over Temp Neon, Buzzer & Reset Button

If your incubator overheats, the buzzer will sound and the over temperature neon will illuminate. **The heaters will automatically disconnect to prevent further overheating.** Once the chamber temperature has stopped rising, and has fallen below the cut out set point, the heaters will re-connect automatically. The buzzer & neon will remain on until a member of staff has acknowledged the alarm. To reset the buzzer & alarm, press the reset button alongside the illuminated neon.

3.6 Setting An Under Temperature Cut Out Value

If, for example, you would like to operate your incubator at +20.0°C, an under temperature cut out value needs to be programmed into the cut out. We recommend that the cut-out value is set to **2 degrees lower** than your chamber temperature. In the example above, the under temperature cut out set point needs to be set +18°C. To do this follow steps 1 - 3 below:

- a) Press the **[P]** button **twice**. Flashing **AL2** will be displayed plus a value (the **under temperature** cut out value).
- b) Press the **[▲]** button to increase the **under temperature** cut out value or the **[▼]** button to decrease the **under temperature** cut out value.
- c) Once the value of your choice is displayed, press the **[P]** button again and the new value will be stored and retained in permanent memory.


3.7 Under Temp Neon, Buzzer & Reset Button

If your incubator overcools, the buzzer will sound and the under temperature neon will illuminate. **The cooling system will automatically disconnect to prevent further overcooling.** The buzzer & neon will remain on until a member of staff has acknowledged the alarm. To reset the buzzer & alarm, press the reset button alongside the illuminated neon.

- 3.8 Allow the chamber temperature to stabilise before putting your products into the chamber.

4 Humidity Control

A wide range of fixed and relative humidities can be achieved by adding a suitable saturated salt solution to the base of the chamber. The relative humidity of air in contact with saturated salt solutions reaches known equilibrium values at certain temperatures.

Salts Graph

The graph above shows the temperatures and relative humidities for several of the more common saturated salts. The distilled water/saturated salt solution should be topped up regularly and should fully cover the floor AND high humidity element at all times.

4.1 High Humidity Option

A high humidity level of approx. 98%RH is maintained by placing distilled water in the base of the chamber and switching the HIGH HUMIDITY switch ON.

Do not use the HIGH HUMIDITY switch unless the floor of the chamber is fully covered with water. It is essential that the heating element on the chamber floor is covered at all times. This must be checked daily.

4.2 Allow the chamber humidity to stabilise before putting your products into the chamber.

5 Maintenance

No routine maintenance is necessary but regular cleaning of the inner chamber is essential.

5.1 Keep ventilation grilles clear and unobstructed.

-
- 5.2 The chamber should be regularly de-frosted as a substantial build-up of ice in the air circulation duct will restrict the airflow and affect the incubators performance. It is recommended that the service ports on the right of the incubator remain sealed with the threaded red bungs supplied to reduce the tendency for humidity in the air to condense out in the chamber.

APPENDIX A

Digital Data Recorder (Temperature & Humidity)

1 Introduction

The Eurotherm nanodac Digital Data Recorder is a simple to use modern version of an analogue paper chart recorder. Data can be viewed on screen or downloaded for viewing on a computer. The digital data recorder has the following features and benefits:

- 3.5" colour TFT screen.
- 50 MB built in non-volatile memory for recording temperature and humidity data.
- Panel mounted USB socket (for downloading recorded data).
- Scrollable colour graphs for **temperature** and **humidity** displayed in real time.
- Data is recorded and stored every few seconds in .csv format
- Recorded data can be imported into Microsoft[®] Excel[®] via a USB flash drive.

-  Recording of data (to built-in memory) starts automatically after power on.
-  The recorder will hold approximately 1 years' worth of data before data is overwritten.
-  Recorded data can be downloaded onto a USB flash drive for analysis (as required).
-  The recorder does not need a USB flash drive inserted to record data. (It only needs one for *downloading* recorded data).

In the event of a power cut, stored data will remain safe in the recorder. Recording will not happen during a power cut.

2 Operation

2.1 Switching On

After plugging the cabinet in, the data recorder will switch on and show a logo while starting up. (A fuel gauge below the logo shows the loading progress).

Please wait for the recorder to start up

The recorder has been factory configured to record data automatically, once it has fully started up. No button pressing required.

The screen (below) will appear once the recorder is operating normally:

Normal operating screen (°C)



Normal operating screen (%RH)

- i** The screen changes automatically between **temperature** and **humidity** every few seconds during normal operation.

2.2 Normal Operation

The recorder starts recording data automatically. The screen shows a live animated graph (in real time) for both the **temperature** and the **humidity**. The screen changes between **temperature** and **humidity** every few seconds. Large numbers in the top right corner of the screen show the actual **temperature** or actual **humidity** so it can be verified easily.


The graph updates itself every few seconds and draws a history line. Historic data from this time line can be viewed from the main menu.

2.3 Control Panel Buttons


It is not necessary to enter the menu or change any settings for normal operation. The recorder has been pre-configured to work optimally without any adjustment or input from the user.

The buttons on the recorders control panel perform these functions:


2.4 Screen Icons

The lower part of the screen contains the status bar, which shows the current date / time. The following icons may also be visible:

A flashing yellow warning triangle means there is no USB flash drive connected.
A white envelope is a system message (used by engineers only).

A green and black R icon means data is being recorded and stored in internal memory (temperature and humidity).

2.5 Viewing Graph Data On Screen

Historical graph data can be viewed on-screen without affecting the recording process.
To view historical graph data:

- 1) Press the  button.


The following screen will appear:

- 2) Press the  button and select **History**.


- 3) Press the  button and a grey background graph screen will be visible. (This is historic graph data that has been recorded for **temperature** and **humidity**).

- 4) Press the  button to scroll BACK in time.

- 5) Press the  button to scroll FORWARD in time.
- 6) Press the  button to change between **temperature** and **humidity** screens.
- 7) When you have finished viewing the historic graph data, press the  button.

The following screen will appear:

- 8) Press the  button and change **No** to **Yes** with the  or  button.
- 9) Press the  button

The following screen will appear:

- 10) Select **Home**.
- 11) Press the  button to exit the Menu.

 The normal operation screen (white background) will be displayed. The screen changes automatically between **temperature** and **humidity** every few seconds.

2.6 USB Socket

A panel mounted USB socket is included next to the recorder. The recorder does not need a USB flash drive inserted to record data - the data is recorded and stored inside the recorder as standard.

The USB socket is only needed when data download to a USB flash drive is required.

A removable screw cap prevents dirt / dust / moisture from entering the USB socket when not in use. Remove the screw cap by turning the cap anti-clockwise to expose the USB socket.

Unscrew cap

Cap removed

Inserting a USB flash drive

Replace the screw cap when you have finished downloading.

2.7 USB Flash Drive Compatibility


The recorder will accept any brand of USB flash drive **up to 8GB** in capacity.

Do not use a USB flash drive larger than 8GB in capacity.

- 1GB ✓
- 2GB ✓
- 4GB ✓
- 8GB ✓
- 16GB ✗
- 32GB ✗
- 64GB ✗
- 128GB ✗
- 256GB ✗

2.8 Downloading Graph Data

The recorder holds approximately 1 years' worth of data internally before data is overwritten. The recorded data is stored in .csv format.

Recorded data can be downloaded onto a USB flash drive when required. This can be done daily, weekly, or monthly.

After the downloaded data has been analysed in a computer, the USB flash drive does not need inserting back into the recorder.

The recorder will continue recording and storing data internally after the USB flash drive has been removed.

When the USB flash drive is next inserted back into the recorder, it will synchronise and update the existing data on that USB flash drive with any new recorded data which has been captured since the last download.

No data acquisition software required for data analysis.

When the recorded data is needed for analysis on a computer, follow this procedure:

Insert a USB flash drive into the USB socket (max. capacity 8GB).

Wait 15 minutes (for the recorded data to synchronise) before removing the USB flash drive.

If it has been a month or longer since you last downloaded recorded data, you may have to wait up to 30 minutes for the data to synchronise.

Remove the USB flash drive from the USB socket.

Insert the USB flash drive into a computer running Microsoft® Excel® (or equivalent).


2.9 Viewing the Downloaded Graph Data

On your computer, navigate to where the downloaded data has been stored on the USB flash drive.

The normal data storage location will be on the flash drive in a folder called “history”. (There may be other sub folders within the “history” folder).

The recorded data will be stored in .csv format in a sub-folder called “CSV”.

An example is shown below:

Open the .csv data in Microsoft® Excel® and display the data in a graph of your own choice.

APPENDIX B

Programming an Offset Parameter into a JUMO dTRON 316 temperature controller

If the temperature displayed by your replacement Jumo 316 controller does not match the actual chamber temperature, an offset parameter may have to be programmed into the controller. To do this, follow the steps below:

1. When the controller is in its normal mode, press the **PGM** button once.
2. **OPr** should be displayed.
3. Press the **_** button twice. **ConF** should appear.
4. When **ConF** is displayed, press the **PGM** button once. **iNP** should appear.
5. Press the **PGM** button twice. **SEnS** should be displayed.
6. Press the **_** button twice. **OFFS°C** should be displayed.
7. To adjust the offset value, press the **PGM** button (**OFFS** flashes), then use the **_** button to increase the OFFS value, or the **_** button to decrease the OFFS value. (If, for example, the temperature display was originally reading 2.0° too high, you need to put a negative offset value in of -2.0. If the temperature display was originally reading 1.0° too low, you need to put a positive offset value in of +1.0).
8. When the desired offset has been entered, **wait a few seconds** and the display will blink to confirm the new offset parameter has been stored.
9. Press the **EXIT** button four times.
10. The current chamber temperature will then be displayed which will include the additional offset value that you have just programmed in.
11. If at any stage during this procedure, the buttons are not pressed for 10 seconds or more, the controller will automatically return to its normal mode and you will have to restart from step 1.

APPENDIX C

TLK38 OFFSET Instructions

1. If the temperature displayed by the TLK38 does not match the calibrated test equipment, the TLK38 can be adjusted by means of an **offset parameter**.
2. Insert a calibrated probe (from your independent test equipment) into the incubator chamber via the access port.
3. Leave the incubator to reach its set point and stabilise.
4. Press and hold the P button for 5 seconds. **ConF** will be displayed.
5. Press the **P** button. **0** will be displayed.
6. Use the **▲** button to make the display read **38.1**
7. Press the **P** button. **SP** will be displayed.
8. Press the **▼** button once only. **1nP** will be displayed.
9. Press the **P** button. **SEnS** will be displayed.
10. Press the **▼** button several times until **OFFSt** is displayed.
11. Press the **P** button.
12. Change the value of **OFFSt** from **0** (factory default) to whatever value is needed to correct the misreading display. (A positive or negative value can be entered).
13. Press the **P** button.
14. Hold the ▼ button for 1 second. **1nP** will be displayed.
15. Press the **▼** button. **OUT** will be displayed.
16. Press the **▼** button. **AL1** will be displayed.
17. Press the **▼** button. **AL2** will be displayed.
18. Press the **▼** button. **LbA** will be displayed.
19. Press the **▼** button. **rEG** will be displayed.
20. Hold the ▼ button for 2 seconds.
21. **The correct temperature should now be showing on the TLK38.**
22. POWER OFF / POWER ON to store the settings and reboot the incubator.

DO NOT CHANGE ANY OTHER SETTING IN THE TLK38 CONTROLLER

DIRECTIVE ON WASTE ELECTRICAL & ELECTRONIC EQUIPMENT (WEEE)

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