



## Pyknometer Complete with Non-Corrodable Cone and Rubber Seal 1 kg Capacity

Code: 24-2885

Product Group: Particle Density: Pyknometer  
Methods, Pycnometer Top and Jar

Glass jar complete with non-corrodable cone.

### Standards

ASTM C128, AASHTO T84

### Accessories



#### 100MI Density Bottle with Perforated Stopper.

Code: 24-2950



#### 50MI Density Bottle with Perforated Stopper.

Code: 24-2900



#### Density Bottle 25MI.

Code: 24-2890



## Pyknometer Gasket

Code: [24-2885/11](#)



## Pyknometer Jar Only

Code: [24-2885/10](#)



## Pyknometer Top

Code: [24-2885/12](#)

## Alternatives



## Sand Absorption Cone

Code: [42-1700](#)