



## Set of Weights 50 kg Slotted.

Code: [26-2132](#)

Product Group: [Application of Normal Stress](#)

Normal stress is applied to the specimen through the lever loading device and the appropriate slotted weights

## Accessories




### Slotted Surcharge Weight 0.25 kg

Code: [26-2147](#)



### Slotted Surcharge Weight 0.5 kg.

Code: [26-2145](#)



### Slotted Surcharge Weight 1 kg.

Code: [26-2143](#)



### Slotted Surcharge Weight 10 kg.

Code: [26-2137](#)



## Slotted Surcharge Weight 2 kg.

Code: [26-2141](#)

---



## Slotted Surcharge Weight 5 kg.

Code: [26-2139](#)

---

### Alternatives



## Slotted Surcharge Weight 0.25 kg

Code: [26-2147](#)

---



## Slotted Surcharge Weight 0.5 kg.

Code: [26-2145](#)

---



## Slotted Surcharge Weight 1 kg.

Code: [26-2143](#)

---



## Slotted Surcharge Weight 10 kg.

Code: [26-2137](#)


---



## Slotted Surcharge Weight 2 kg.

Code: [26-2141](#)

---



## Slotted Surcharge Weight 5 kg.

Code: [26-2139](#)